Sooper Products

INTRODUCES

SooperMask

Intelligent Protective Gear

**My mask saves you,
your mask saves me! #SooperMasks**

Innovative Hygiene Technology for Maximum Protection.

One virus has completely changed the way we live. Making health and hygiene, everyone's top priority. In today's times, people across the globe are anxious to stay safe from the dreaded Covid-19 virus, and are in search of everyday solutions.

Sooper Products offers a range of solutions that are created using innovative hygiene technology to safeguard people. Whether you want to breathe easy from behind a mask when you step out or disinfect surfaces and objects with UV light, we have you covered. With **Sooper Products** on your side, you can rest assured that your health risk is significantly reduced. And so is your stress.

The **SooperMasks** have 3 layers but it gives 5 levels of protection against SARS-CoV-2, the virus which causes Covid-19. Specially treated textiles along with the filter channels eliminate 99.9% of the virus. **Livinguard Technology is applied to the fabric covering the filter media.**

The **SooperLamp** UV box uses the power of UV light to keep surfaces and objects virus-free for you.

As we go about our daily lives, knowing fully well that the fight with virus is going to be a long drawn one, the best we can do is stay highly protected always with **Sooper Products**.

**Destroys
Coronavirus***

SooperMask

Intelligent Protective Gear

Medical Face Mask Type II
Bacteria Filtration Efficacy (BFE) >98%

Reinventing protection by **destroying SARS-CoV-2**

Livinguard Technology is applied to the fabric covering the filter media.

Reusable

washable up to 30 times
and reusable for 6
months.

Sustainable

replaces 210
conventional masks,
does not leach.

Protection

fabric destroys bacteria, fungi
and viruses and prevents
them from spreading.

Style

Available in Vantablack,
Bombay Blue, Forest
Green, Cosmic Red

Comfort

Silky smooth fabric,
odourless, adjustable
ear loops for custom
fit, high breathability

Safety

High performance
filtration & free of
toxic metals

Engineered and designed in Switzerland

Standards of Conformity

SooperMask
Intelligent Protective Gear

Proven Antiviral Properties:

| Type of Virus | Test Facility | Viral Inactivation % |
|------------------------|--|----------------------|
| SARS-CoV-2 | Freie Universität Berlin, Germany | 99.9% |
| Human Coronavirus 229E | University of Arizona, Tucson, USA | 99.5% |
| Influenza A | Federal Office for Civil Protection, FOCP, Switzerland | > 99.96% |
| Yellow Fever | Federal Office for Civil Protection, FOCP, Switzerland | > 99.97% |
| Bacteriophage MS2 | BTS, Mumbai, India | > 99.99% |

Antibacterial against a wide range of bacteria,
for example:

Escherichia coli, *Staphylococcus aureus*, *Candida albicans*, Multidrug-resistant *Enterococcus faecium* (VRE - ATCC 51559), Methicillin-resistant *Staphylococcus aureus* (MRSA - ATCC 33591), *Mycobacterium tuberculosis* (ATCC 25177), *Klebsiella pneumoniae*, *Clostridium difficile* (*C. diff*), *Mycobacterium terrae*, *Mycobacterium avium*, *Salmonella Typhi*, *Pseudomonas aeruginosa*

Why our mask?

Infinitely better and >25x cheaper masks

vs.

Conventional N95

Destroys viruses* and bacteria upon contact without leaching poisonous chemicals / metals used in traditional antimicrobial treatments

Virus & bacteria impact

Viruses survive and bacteria multiply on the mask, thereby increasing infection risk

210 uses based on daily extended use x 30 weekly washes

Usable life

Single use, works for < 8 hours

< \$0.12 per use based on \$25 per mask

Consumer's cost / use

Normally \$1 per mask and per use, *not accounting for current prices*

N95, can make N99 at slightly additional cost

Filtration performance

N95

< 0.5% waste.
No leaching.

Environmental impact

> 200x non-biodegradable N95 filters to be disposed.

Virus cannot mutate or develop resistance.

Poisonous leaching if treated with conventional antimicrobials.

* See Appendix for virus types and test results

Livinguard patented solution for respirators

Layer 2 (N95) blocks particulate matter


- Nonwoven filtration media blocks most particles in the air and therefore most of the germs contained in dust

Livinguard-charged layers (1 and 3) destroy microbes

- Microbes on dust particles filtered by layers 1 and 3 are partially killed (in contact with textile surface)
- Microbes contained in droplets filtered / absorbed by layer 1 and 3 are killed very effectively
- Microbes accumulating in layer 2 and transferred through respiration to layers 1 and 3 are killed upon contact with the textile surface
- Layers 1 and 3 remain disinfected and show very low contamination levels

Results

- Strongly reduces the transfer of pathogens from infected persons to the environment and from the environment to non-infected persons
- As the outside layers (1 and 3) are permanent disinfectants, **cross-contamination of the person handling the mask is largely reduced**



Why our mask?

Proven efficacy against other Viruses

NSF International (USA) test report J-00293323
BTS (India ISO & NABL accredited lab) test report 14386/1-3

FR coliphage: 99.999% reduction (5 log)


MS2 bacteriophage: > 99.999% reduction (> 5 log)
from NSF and > 99.99% (> 4 log) from BTS



Proven efficacy against other Viruses

Testing by Swiss Federal Lab ABC Spiez

Yellow Fever Virus:
> 99.97% reduction (> 3.6 log)

Influenza A Virus:
>99.96% reduction (> 3.5 log)

How does Livinguard destroy Viruses?

- 1 cm² of Livinguard treated textiles has an array of 24 billion immobilized cationic charges of different amplitudes
- These positive charges interact with and extract negatively charged compounds from the cell membrane (bacteria) or capsid (viruses), leading to disruption of the integrity and inactivation of the organism, without leaching


Image reference of an exploded Virus using cationic charges:
[B.B. Hsu, S.Y. Wong, P.T. Hammond, J. Chen, A.M. Klibanov, PNAS January 4, 2011 108 \(1\) 61-66](#)

Certifications

China National Accreditation Service for Conformity Assessment INSPECTION BODY ACCREDITATION CERTIFICATE

(Registration No. CNAS IB0024)

Guangdong Testing Institute of Product Quality Supervision

No.10, Science Avenue, Science City, Huangpu District, Guangzhou,
Guangdong, China.

is accredited in accordance with ISO/IEC 17020:2012 General Criteria for the Operation of Various Types of Bodies Performing Inspection (CNAS-CI01 Accreditation Criteria for the competence of Inspection Bodies) as Type A inspection body for the competence to undertake inspection services as described in the schedule attached to this certificate.

The scope of accreditation is detailed in the attached schedule bearing the same registration number as above. The schedule forms an integral part of this certificate.

Date of Issue: 2018-01-03

Date of Expiry: 2024-01-15

Date of Initial Accreditation: 2004-09-13

Signed on behalf of China National Accreditation Service for Conformity Assessment

China National Accreditation Service for Conformity Assessment (CNAS) is authorized by Certification and Accreditation Administration of the People's Republic of China (CNCA) to operate the national accreditation schemes for conformity assessment. CNAS is a signatory of the International Laboratory Accreditation Cooperation Mutual Recognition Arrangement (ILAC-MRA) and the Asia Pacific Laboratory Accreditation Cooperation Mutual Recognition Arrangement (APLAC-MRA). The validity of the certificate can be checked on CNAS website at <http://www.cnas.org.cn/english/findanaccreditedbody/index.shtml>

SooperMask

Intelligent Protective Gear

EC declaration of conformity

This is to declare that:

Livinguard AG
Bahnhofstrasse 12
6300 Zug
Switzerland

Product Name:

Livinguard Antibacterial & Antiviral N95/N99 Reusable Face Mask with/without valves

Product model:

FM-95-30-WV-BV: Face mask / N95 / 30 washes / with valve / Antibacterial & Antiviral
FM-95-30-NV-BV: Face mask / N95 / 30 washes / no valve / Antibacterial & Antiviral
FM-99-30-WV-BV: Face mask / N99 / 30 washes / with valve / Antibacterial & Antiviral

Livinguard respiratory face masks Livinguard respiratory face mask can filter >95% droplets, airborne and particles (>0.3 microns). The face masks are hygienic, antibacterial, antiviral, reusable, and washable.

The manufacturer bears the sole responsibility for issuing this declaration of conformity.

Livinguard face masks comply with the Regulation (EU) 2016/425 on personal protective equipment for respiratory protection.

Harmonized standards applied in whole or in part: EN 149:2001+A1:2009 (Reusable FFP2).

EC Declaration of Conformity: As an authorized representative, I declare that the above mentioned products meet the requirements of the applicable directives.

Signed for and on behalf of:

Mr. Sanjeev Swamy CXO

p.a. Wouter Prok, EVA R&D

Photo verification of the above mentioned object:

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The observation of perspiration extract showed no discoloration. In the investigation on odor a slight sweet smell was detected.

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The pH-value in the perspiration extract was still 5.5.

| Rel. Proteingehalt / Rel. protein content OD [570nm] | X | ± | s | Wachstumshemmung / growth inhibition in % |
|--|--------|---|--------|---|
| Leerwert / blank: | 0,1894 | | | |
| Positivkontrolle / positive control: | 0,3719 | ± | 0,0079 | 88 |
| Negativkontrolle / negative control: | 1,9360 | ± | 0,0832 | 0 |
| Lösungsmittelkontrolle / solvent control: | 1,6762 | ± | 0,0041 | 0 |
| Prüfgegenstand / test sample: | | | | Limit |
| 33,30% | 1,4613 | ± | 0,0625 | 14 |
| 22,20% | 1,5517 | ± | 0,0361 | 8 |
| 14,80% | 1,6224 | ± | 0,0586 | 4 |
| 9,90% | 1,6163 | ± | 0,0206 | 4 |

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X is the average of multiple measurements
s is the standard deviation
Test performance: hec, nka

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Under the mentioned test conditions, the perspiration extract of the test material showed a growth inhibition of 14 % in the cytotoxicity test.

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A dose-dependent growth inhibition of L929 cells could be observed, which derived from the test material, but the significant limit of 30 % was not exceeded.

Australian Government
Department of Health
Therapeutic Goods Administration

Australian Register of Therapeutic Goods Certificate

Issued to

M4 Medical Pty Ltd

for approval to supply

M4 Medical Pty Ltd - Mask, surgical, reusable

| | |
|------------------|--|
| ARTG Identifier | 335554 |
| ARTG Start Date | 1/05/2020 |
| Product Category | Medical Device Included Class 1 |
| GMDN | 37713 |
| GMDN Term | Mask, surgical, reusable |
| Intended Purpose | An N95 reusable mask placed over the nose and mouth by medical personnel to prevent the transmission of airborne organisms while surgery is being performed. |

| Manufacturer Details | Address | Certificate number(s) |
|-------------------------|--|-----------------------|
| FINE HYGIENIC PAPER FZE | PLOT NO EWRA 56 JEBAL ALAI FREE ZONE PO BOX 16927 , DUBAI , United Arab Emirates | |

ARTG Standard Conditions

The above Medical Device Included Class 1 has been entered on the Register subject to the following conditions:

- The inclusion of the kind of device in the ARTG is subject to compliance with all conditions placed or imposed on the ARTG entry. Refer Part 4-5, Division 2 (Conditions) of the Therapeutic Goods Act 1989 and Part 5, Division 5.2 (Conditions) of the Therapeutic Goods (Medical Devices) Regulations 2002 for relevant information.
- Breaching conditions of the inclusion related to the device of the kind may lead to suspension or cancellation of the ARTG entry; may be a criminal offence; and civil penalties may apply.

Products Covered by This Entry

1. Mask, surgical, reusable

Product Specific Conditions

No specific conditions have been recorded against this entry.

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According to DIN EN ISO 10993-5 a growth inhibition of more than 30 % in comparison to the solvent control is considered as a clear cell-toxic effect.

BEURTEILUNG

ASSESSMENT

Probe / Sample 18.8.5.0138-3

Extrakt von Reusable cotton hand towel roll scoured treated (1 x washed by Hohenstein Institute)

Extract of Reusable cotton hand towel roll scoured treated (1 x washed by Hohenstein Institute)

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The detected growth inhibition of less than 30 % shows that no cytotoxic substances are released.

Schloss Hohenstein, 23. Februar 2018

Deputy Director
Life Science & Care

Anja Gerhardt

Dr. Anja Gerhardt

Product Manager Biomedicine
Life Science & Care

Anne Schopfer

Dipl.-Ing. (FH) Anne Schopfer

"The translation was carried out to the best of a non-native speaker's knowledge. Liability cannot be taken."

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RESEARCH, DISCOVERY & INNOVATION

Water & Energy Sustainable Technology Center

Study Conclusions

The FM-7 Test Fabric was the most efficacious against Human Coronavirus 229E, achieving a reduction of 99.5% (2.30 log₁₀) given two hours of exposure relative to the Control Fabric (Table 1). Test Fabric FM-26 reduced levels of the test virus by 98.4%, while reductions were minimal for Test Fabric FM-14 at 13.8%.

With regard to cytotoxicity effects on the MRC-5 cell line, the level of toxicity observed for test fabrics FM-7 and FM-26 was similar to that of the Control Fabric (Table 2). Toxicity levels on MRC-5 host cell monolayers were slightly greater for the FM-14 test fabric. However, the sensitivity of MRC-5's, an internal human lung fibroblast cell line, should not necessarily be extrapolated assume any effects on the protective and resilient cells comprising the external cell layers of human skin.

Neutralization was validated for each of the Control and Test Fabrics beyond the level of cytotoxicity, with < 0.5 log₁₀ difference in viral titer observed between the Control and Test fabrics (Table 3).

References

- Sturman, L.S., Holmes, K.V. and Behnke, J., 1980. Isolation of coronavirus envelope glycoproteins and interaction with the viral nucleocapsid. *Journal of virology*, 33(1), pp.449-462.
- Chen, D., Luo, W., Hoffman, J., Huang, L., Sandefur, S., Hall, T., Murphy, M. and O'Donnell, S., 2019. Insights into Virus Inactivation by Polysorbate 80 (PS80) in the Absence of Solvent. *Biotechnology Progress*.
- Colavita, F., Quartu, S., Lalle, E., Bordini, L., Lapa, D., Meschi, S., Vulcano, A., Toffoletti, A., Bordini, E., Paglia, M.G. and Di Caro, A., 2017. Evaluation of the inactivation effect of Triton X-100 on Ebola virus infectivity. *Journal of Clinical Virology*, 86, pp.27-30.

Biocompatibility of Mask Material
Prepared For: Livinguard, Inc.
22 April, 2020

Biocompatibility of Mask Material
Prepared For: Livinguard, Inc.
22 April, 2020

Livinguard Self-Disinfecting Textile for Emergency Use in N95 Facemasks: Statement of Biocompatibility

Background

COVID-19 is an infectious disease caused by the SARS-COV-2 virus that has been declared a global pandemic by the World Health Organization. The spread of SARS-COV-2 has resulted in immense stress on global healthcare systems and produced a shortage of lifesaving medical devices and personal protective equipment used by medical professionals, including face masks. Livinguard, Inc. manufactures a range of self-disinfecting materials including textiles that are fabricated through the integration of disinfectants directly into the material. Livinguard wishes to use their self-disinfecting textile in the fabrication of an N95 facemask. This use would classify the device as an N95 Respirator with Antimicrobial/Antiviral Agent (product code ONT) regulated by 21CFR878.4040 which specifies that "the user contacting components of the device must be demonstrated to be biocompatible." The FDA has raised a specific question on the point of biocompatibility of the treated textile when used as a mask.

The textile in question is treated with polyhexamethylene biguanide (PHMB, Manufactured by Arch UK Biocides as Vantocil IB). This active agent is bonded to the textile fibers, with effectiveness that persists after multiple washes (data on file with Livinguard). This document summarizes what is known about the biocompatibility of the Livinguard self-disinfecting textile in question, and evaluates the appropriateness of its use in N95 masks from the perspective of biocompatibility.

Summary of Biocompatibility Testing Completed on MED610

A suite of biocompatibility tests have been conducted on a worst-case version of the textile that contains PHMB among other ingredients at a concentration 250% above normal.

Table 1. Biocompatibility Tests Performed on Livinguard's Self-Disinfecting Textile

| Biological Endpoint | Testing Standard | Result |
|-----------------------|-------------------------------------|--------|
| Sensitization | EPA P328 Buehler Method | Pass |
| Irritation | EPA P326 Primary Skin Irritation | Pass |
| Acute Dermal Toxicity | EPA P322 Acute Dermal Toxicity | Pass |

Gap Assessment, Testing Conducted, and Regulatory Expectations for Mask Applications

Testing was conducted on the textile per EPA protocols assuming household use, not necessarily as a mask. For masks that are marketed as a medical device, a biocompatibility evaluation should be conducted per ISO 10993-1. According to this standard, N95 masks would be categorized as surface devices with a cumulative contact to intact skin possibly greater than 30 days. The biological endpoints for evaluation of this contact type are cytotoxicity, sensitization, and irritation.

Based on the passing results of sensitization and irritation per EPA standards, it is likely that the material would pass these tests when tested per ISO 10993-1, and the risks associated with these endpoints are acceptable. Cytotoxicity testing was not conducted. This test would not be a good indicator of biocompatibility for the textile in question, as it is designed to be anti-microbial and is expected to be cytotoxic. Cytotoxicity testing should be considered on the same textile not treated with the active agent. Nevertheless, the absence of these test results should not prevent the material from being used in PPE during the COVID-19 pandemic.

FDA Feedback Regarding Use

In preliminary conversations, the FDA indicated to Livinguard a concern that the charged fibers of the textile could act as a cationic surfactant that may damage respiratory tissue.

Vantocil IB is an aqueous solution of PHMB, which is applied to the textile and then fully dried. PHMB is a polymer, typically with 10-40 repeating units and a molecular weight in excess of 700 g/mol. Therefore, it is expected to be non-volatile and has a reported vapour pressure of 4.11x10⁻⁷ Pa at 25 °C and not able to vaporize for inhalation at any appreciable amounts (see CLH Monograph on PHMB). Because the active agent on the treated textiles is permanently applied, persisting over several washes, and non-volatile they are not available to be inhaled and interact with respiratory tissues.

Potential Risks Associated with VOCs and Particulates: While VOCs released from the textile were not explicitly measured, the non-volatile nature of active agents and their persistence on the textile mitigate this concern. Therefore, any risks potentially associated with VOCs and particulates from the textile are considered acceptable in emergency situations associated with the COVID-19 pandemic.

Limits and Recommendations

The risks related to use of Livinguard's self-disinfecting textile in masks is considered minimal. It is recommended that representative samples of material are fabricated into masks and tested per the most current FDA recommendations to explicitly demonstrate compliance to the standards and to completely mitigate risk. It is however acceptable to proceed with use in applications related to the COVID-19 pandemic ahead of the conclusion of this testing.

Conclusion

The minimal biocompatibility risks associated with use of Livinguard self-disinfecting textile in masks are acceptable in emergency situations associated with the COVID-19 pandemic.

Document Prepared By:

Matthew R. Jorgensen

2020.05.02 10:14:59
-06'00'

Matthew R. Jorgensen, PhD, DABT
Chemistry and Materials Scientist
Nelson Laboratories, LLC
P: 801-290-7794
E: mjorgensen@nelsonlabs.com
NELSON LABORATORIES, LLC
A Sotera Health Company

Important Information: This document is not equivalent to regulatory approval by any regulatory body. It is a third-party expert evaluation of the potential risks of the material when used described and in accordance with Emergency Use Authorization (EUA) being issued in response to concerns relating to insufficient supply and availability of FDA-cleared PPE for use in healthcare settings to treat patients during the Coronavirus Disease 2019 (COVID-19) pandemic. For more information visit US FDA at: <https://www.fda.gov/media/136423/download>.

Certifications

Supplementary Information about Chinese Mask Test Report (FZ2011002)

27.04.2020, Zug

Livinguard Respiratory Face Mask was tested by Chinese lab GUANGDONG TESTING INSTITUTE OF PRODUCT QUALITY SUPERVISION, the test results show Livinguard Respiratory Face Mask is compliant with Chinese Standard GB 2626 2006 KN95 (Respiratory protective equipment. Non-powered air-purifying particle respirator), based on the comparison below, it is reasonable to consider China KN95 as "equivalent" to US NIOSH N95 and European FFP2 respirators.

| Certification/Class (Standard) | N95 (NIOSH-42C FR84) | FFP2 (EN 149-2001) | KN95 (GB2626-2006) | P2 (AS/NZ 1716:2012) | Korea 1st Class (KMOEL - 2017-64) | DS (Japan JMHLWNotification 214, 2018) |
|--|----------------------|--------------------------------|--------------------------------|---|-----------------------------------|---|
| Filter performance – (must be ≥ X% efficient) | ≥ 95% | ≥ 94% | ≥ 95% | ≥ 94% | ≥ 94% | ≥ 95% |
| Test agent | NaCl | NaCl and paraffin oil | NaCl | NaCl | NaCl and paraffin oil | NaCl |
| Flow rate | 85 L/min | 95 L/min | 85 L/min | 95 L/min | 95 L/min | 85 L/min |
| Total inward leakage (TIL)* – tested on human subjects each performing exercises | N/A | ≤ 8% leakage (arithmetic mean) | ≤ 8% leakage (arithmetic mean) | ≤ 8% leakage (individual and arithmetic mean) | ≤ 8% leakage (arithmetic mean) | Inward Leakage measured and included in User Instructions |

Livinguard AG
Bahnhofstrasse 12, 6300 Zug, Switzerland

| Certification/Class (Standard) | N95 (NIOSH-42C FR84) | FFP2 (EN 149-2001) | KN95 (GB2626-2006) | P2 (AS/NZ 1716:2012) | Korea 1st Class (KMOEL - 2017-64) | DS (Japan JMHLWNotification 214, 2018) |
|---|-----------------------|--|-----------------------------------|---|---|---|
| Inhalation resistance – max pressure drop | ≤ 343 Pa | ≤ 70 Pa (at 30 L/min) ≤ 240 Pa (at 95 L/min) ≤ 500 Pa (clogging) | ≤ 350 Pa | ≤ 70 Pa (at 30 L/min) ≤ 240 Pa (at 95 L/min) | ≤ 70 Pa (at 30 L/min) ≤ 240 Pa (at 95 L/min) | ≤ 70 Pa (w/valve) ≤ 50 Pa (no valve) |
| Flow rate | 85 L/min | Varied – see above | 85 L/min | Varied – see above | Varied – see above | 40 L/min |
| Exhalation resistance – max pressure drop | ≤ 245 Pa | ≤ 300 Pa | ≤ 250 Pa | ≤ 120 Pa | ≤ 300 Pa | ≤ 70 Pa (w/valve) ≤ 50 Pa (no valve) |
| Flow rate | 85 L/min | 160 L/min | 85 L/min | 85 L/min | 160 L/min | 40 L/min |
| Exhalation valve leakage requirement | Leak rate ≤ 30 mL/min | N/A | Depressurization to 0 Pa ≥ 20 sec | Leak rate ≤ 30 mL/min | visual inspection after 300 L /min for 30 sec | Depressurization to 0 Pa ≥ 15 sec |
| Force applied | -245 Pa | N/A | -1180 Pa | -250 Pa | N/A | -1,470 Pa |
| CO2 clearance requirement | N/A | ≤ 1% | ≤ 1% | ≤ 1% | ≤ 1% | ≤ 1% |

*Japan JMHLW-Notification 214 requires an Inward Leakage test rather than a TIL test.

Livinguard AG
Bahnhofstrasse 12, 6300 Zug, Switzerland

Berlin, 16.06.2020

Report-Antiviral Activity of textile products from Livinguard

Tested by: Prof. Uwe Rösler and Dr. Anika Friese

Freie Universitaet Berlin, Institute for Animal Hygiene and Environmental Health

Tested textile product

- reference textile: Baumwolle (Co) unbehandelt 11.04.2020, lab number: C1
- antiviral textile: HEFM47.CO.01, 13.05.2020, LIVINGUARD, lab number: C2-neu

Method

- based on ISO 18184 (First edition 2014-09-01)
 - washing of all tested textiles with deionized water 10 times at 40°C, drying
 - cutting pieces of approximately 20x20mm and making up a mass of 0,4g with several pieces
 - sterilization at 121°C for 15 min, drying
 - before starting the test **conditioning the textiles overnight** in a humid environment (incubator 37°C)
 - controls: Verification of cytotoxic effect and cell sensitivity to virus/inactivation of antiviral activity
 - test: inoculation of 0,4g textile **with 1ml virus suspension** (at least 10⁷ TCID₅₀/ml) washing out using 19 ml **cell culture medium** auf (MEM Eagle EBSS + 10% FKS + 1% nicht-essentielle AS, +1 % Penicillin (10.000 U/ml) / Streptomycin (10 mg/ml)) titration in 96-well plate and titer calculation with method of **Spearman and Karber in TCID₅₀/ml**
 - comparison of titer from untreated textile (reference) and treated textile for the specific time points

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Evaluation of the conformity of corona sars-cov-2 pandemic respiratory protection (CPA) according to the testing principle for corona sars-cov-2 pandemic respiratory protection masks revision 1

| | |
|-------------------------------------|--|
| Berichtsnummer <i>Report number</i> | 8003019204-1 |
| Prüfgegenstand <i>Test subject</i> | Partikelfiltrierende Halbmaske ohne Ventil |
| Modell <i>Type</i> | High-Performance Virus-Eliminating Mask |
| Importeur <i>Importer</i> | Medivid GmbH Billwerder Neuer Deich 70 20539 Hamburg |
| Hersteller <i>Manufacturer</i> | Anhui Tiankang Medical Technology Co. Ltd. 228 Weyi Road, Economic Development Zone, Tianchang City, 239300 Anhui, China |

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SooperMask

Intelligent Protective Gear

Suitable for daily personal use.
Medical Face Mask Type II

- Conforms to EN14683 Medical Mask Type II (Bacteria Filtration Efficacy (BFE) >98%)
- CE Medical device class I
- Choice of 4 solid colours

Destroys
99.9%
Coronavirus

Also available for kids:

SooperMask

Intelligent Protective Gear JUNIOR

Contact us for more details.

SOOPER PRODUCTS

hiren@sooperproducts.com

Hiren Kakad +919820331000

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TÜV NORD CERT GmbH
TÜV NORD CERT GmbH, Langemarkstraße 20, 45141 Essen